

CareFIT™ P2 / N95



3-PANEL MEDICAL GRADE RESPIRATOR

MSK-005



- Australian-made, Made in Geelong
- Australian-certified P2 Respirator AS/NZS 1716:2012
- Tested to N95 NIOSH standard 42 CFR Part 84 subpart K



FEATURES



THREE PANEL DESIGN

Proprietary three-panel design delivers structural integrity, an effective facial fit and optimal comfort.




SPLASH BARRIER

Provides the highest-recognised level of splash protection. Synthetic blood penetration resistance: 160 mmHg (Level 3 Fluid Resistance).




ULTIMATE PARTICULATE PROTECTION

The electrostatically charged fibres in the filter media creates a tortuous path which captures $\geq 95\%$ of all particles, including the most difficult at 0.3 microns when the respirator is securely fitted.



HIGH INDIVIDUAL FIT

The adjustable foam nose pad provides a secure seal that is adaptable to the contours of the face to support high individual fit. An effective design for many 'hard to fit' individuals.



OPTIMAL USER COMFORT

The inner lining layer is soft and breathable, imparting a low breathing resistance to maximise wearer comfort.



AUSTRALIAN MADE

TGA REGISTERED ARTG No.: 340518

Manufactured by Care Essentials Pty. Ltd.

103 Mornington Street,
North Geelong, VIC, Australia 3215
www.careessentials.com.au

Customer Service

T +61 3 5277 1455
E queries@careessentials.com.au