

by **Desalys** 

5

STERIPLUS<sup>TM</sup> 80

ADVANCED BIOMEDICAL WASTE TREATMENT SOLUTIONS



## **STERIPLUS™80**

INTEGRATED SYSTEM WITH PRELIMINARY SHREDDING AND STEAM STERILIZATION OF YOUR BIOHAZARDOUS WASTE

The STERIPLUS<sup>™</sup> systems from Tesalys are the ideal solution to inactivate your biohazardous waste safely on-site. Its integrated shredding system TESASHRED<sup>™</sup> not only reduces the volume and mass of waste, but also ensures waste is fully processed.

] () N

0V

A

**TESALYS** Main address 7, rue du Cassé 31240 Saint-Jean (Toulouse) France

Tel. +33 5 62 10 18 91 info@tesalys.fr

www.tesalys.fr





## STERIPLUS<sup>™</sup>80

integrated system with preliminary shredding and steam sterilization of your biohazardous waste



For consumption with cooling of the drain and / or the Virus + kit, please consult us.

## **STERIPLUS<sup>™</sup>** A COMPLETE RANGE



STERIPLUS<sup>™</sup> 20

2/5 kg/h



STERIPLUS<sup>™</sup> 40 4/10 kg/h



STERIPLUS<sup>™</sup> 80 8/20 kg/h

**USB PORT** 



Built-in

**TESABAG** Autoclavable biohazardous waste bags



**TESABOX** 

**Cardboard containers** 

for biohazardous waste

vithout They are subject to change and are not binding. -Ref. this brochure are for Illustrations: © Clément Debeir ation, photos and illustrations in digeek.fr — hnical informa